CellCor™ EXO CD User Guide

Serum-free chemically defined media for expansion and isolation of hMSC-derived exosome

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Product No.	Product	Volume	Storage	Shelf life
YSP017	CellCor™ EXO CD	500 mL	2-8°C	6 months

Key Requirements

항 목	권 장 사 항
Antibiotics	Gentamicin or Penicillin streptomycin
Detach solution	TrypLE™ Express (Gibco)
Culture ware	Tissue Culture Flask, Plate, Cell Factory (Corning, Nunc or Falcon)
Seeding cell counts	4,000-5,000 cells/cm ²
Media volume	T25 flask (5 mL), T75 flask (15 mL), T175 flask (30-40 mL), Cell Factory (100mL-1L)
Cell harvest confluency	75-85% (Figure 참조)
Exosome isolation confluency	85-90% (Figure 참조)

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Process (at a glance)

Enabling Flexible Research Design

■ Media Change Process

Cell thawing with CellCor™ MSC CD AOF (or commonly used culture media) → Cell expansion → (Rinsing) → Media Change with CellCor™ EXO CD → Exosome enrichment → Conditioned media Collection → Exosome isolation

■ Seeding with CellCor™ EXO CD Process

Cell thawing with CellCor™ MSC CD AOF (or commonly used culture media) → Cell expansion → Cell seeding with CellCor™ EXO CD → Exosome enrichment → Conditioned media Collection → Exosome isolation

■ Thawing with CellCor™ EXO CD Process

Cell Thawing with CellCor™ EXO CD → Cell expansion → Exosome enrichment → Conditioned media Collection → Exosome isolation

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