Product	Catalog #	Grade	Size	Storage	Shelf Life*
CellCor™ SFD MSC	YSP010	GMP	500 mL	Under -20°C	Expiry date on label

* Shelf Life is 12 months from Date of Manufacture.

Product Description

CellCor™ SFD MSC is a Serum-free defined medium for expansion and growth of human mesenchymal stem cells (hMSCs, i.e. derived from adipose tissue, bone marrow, umbilical cord blood; AdMSC, BMSC, UCSC). CellCor™ SFD MSC does not contain serum such as fetal bovine serum (FBS). This product conducts Quality Control tests including sterility, mycoplasma, and endotoxin test.

※ Coating agent may be used if needed.

※ This product does not contain antibiotics.

Intended Use

For Research Use or Further Manufacturing Use Only. Not intended for diagnostic use or as for direct therapeutic human and veterinary use.

Safety Information

Read and follow the instructions on the Material Safety Data Sheet (MSDS). Be sure to wear appropriate protective eyewear, clothes and gloves.

Medium Preparation

1. CellCor™ SFD MSC should be thawed at room temperature for 1 hour, then completely thawed in a 37°C water bath.
2. Thawed CellCor™ SFD MSC is stable at 4°C for 4 weeks.
3. Before using CellCor™ SFD MSC for cell culture, aliquot and warm the medium in a 37°C water bath for 30 minutes.
 - ※ Use immediately after thawing. Avoid additional freeze-thaw cycles.

hMSC Culture Protocol

Cell Type: human mesenchymal stem cells (i.e. AdMSC)

Culture Condition: 37°C in a humidified atmosphere of 5% CO₂

Culture Vessel: T-75 flask (i.e. Corning, Catalog #430641U)

[Recovery]

1. Disinfect the interior of BSC with 70% ethanol, and prepare warmed CellCor™ SFD MSC. Add 9 mL of pre-warmed CellCor™ SFD MSC into a 15 mL tube.
2. Thaw cryovial in a 37°C water bath for 1-2 minutes.
3. Disinfect the surface of the cryovial with 70% ethanol and place it in the BSC.
4. Carefully transfer the thawed cells to the 15 mL tube containing CellCor™ SFD MSC.
5. Centrifuge the tube at 230 x g, 20°C for 3 minutes.
6. Remove the supernatant and suspend cells with CellCor™ SFD MSC. Count cells and seed to flask at density 4,000-5,000 cells/cm² (i.e. AdMSC) with 15 mL of CellCor™ SFD MSC.
7. Incubate at 37°C in a humidified atmosphere of 5% CO₂.

[Subculture]

1. Remove the cultured medium and wash with 10 mL of PBS.
2. Add 2 mL of Accutase (Life Technologies, Catalog #A1110501) or TrypLE Express (Life Technologies, Catalog #12604013) to the flask and incubate cells for 1-3 minutes at room temperature.
 - ※ When using detachment reagents such as Trypsin, appropriate inhibitors such as Trypsin inhibitors are required to inactivate enzymes.
3. Observe the cell under the microscope for cell detachment.
 - ※ Gently tap the flask when cells do not detach.
4. Using CellCor™ SFD MSC, collect cells and transfer to the tube.
5. Centrifuge cells at 230 x g, 20°C for 3 minutes.
6. Remove the supernatant and suspend the cell pellets with CellCor™ SFD MSC. Count cells and seed to flask at density 4,000-5,000 cells/cm² (i.e. AdMSC) with 15 mL of CellCor™ SFD MSC.
7. Incubate at 37°C in a humidified atmosphere of 5% CO₂.
8. Subculture when cells reach up to 75-85% confluence. (Figure 1, usually 3-4 days post seeding)


Figure 1. Cell Confluency (75-85%)

[Cryopreservation]

1. Remove the supernatant and suspend cell pellets with 1 mL of cryopreservation medium (e.g. Cell-banker2, Masbio, Catalog #11891).
2. Immediately place cryovials in a freezing container and place at -80°C for overnight. After 24 hours, place and store cryovials in a LN₂ tank.

[Adaptation from medium to CellCor™ SFD MSC]

1. After culturing with a prior medium for 24 hours, remove half of the medium and replace the same amount with CellCor™ SFD MSC.
 - ※ i.e. In a T-75 flask, Culture with 15 mL of the prior medium for 24 hours. Remove 7.5 mL of prior medium and add 7.5 mL of CellCor™ SFD MSC.
2. When cells reach up to 75-85% confluency, subculture using CellCor™ SFD MSC.

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