# CellCor™ CD MSC (RUO)

Product	Catalog #	Grade	Size	Storage	Shelf Life*
CellCor™ CD MSC	YSP001	RUO	500 mL	Under -20°C	Expiry date on label

\* Shelf Life is 12 months from Date of Manufacture.

## Product Description

CellCor™ CD MSC is a Serum-free chemically defined medium for expansion and growth of human mesenchymal stem cells (hMSCs, i.e. derived from adipose tissue, bone marrow, umbilical cord blood; AdMSC, BMSC, UCSC). CellCor™ CD MSC does not contain animal/human nor plant-derived materials. This product conducts Quality Control tests including sterility, mycoplasma, and endotoxin test.

※ Coating agent may be used if needed.

※ This product does not contain antibiotics.

## Intended Use

For Research Use Only. Not intended for human or animal diagnostic or therapeutic uses.

## Safety Information

Read and follow the instructions on the Material Safety Data Sheet (MSDS). Be sure to wear appropriate protective eyewear, clothes and gloves.

## Medium Preparation

1. CellCor™ CD MSC should be thawed at room temperature for 1 hour, then completely thawed in a 37°C water bath.
2. Thawed CellCor™ CD MSC is stable at 4°C for 4 weeks.
3. Before using CellCor™ CD MSC for cell culture, aliquot and warm the medium in a 37°C water bath for 30 minutes.  
※ Use immediately after thawing. Avoid additional freeze-thaw cycles.

## hMSC Culture Protocol

Cell Type: human mesenchymal stem cells (i.e. AdMSC)

Culture Condition: 37°C in a humidified atmosphere of 5% CO<sub>2</sub>

Culture Vessel: T-75 flask (i.e. Corning, Catalog #430641U)

## [Recovery]

1. Disinfect the interior of BSC with 70% ethanol, and prepare warmed CellCor™ CD MSC. Add 9 mL of pre-warmed CellCor™ CD MSC into a 15 mL tube.
2. Thaw cryovial in a 37°C water bath for 1-2 minutes.
3. Disinfect the surface of the cryovial with 70% ethanol and place it in the BSC.
4. Carefully transfer the thawed cells to the 15 mL tube containing CellCor™ CD MSC.
5. Centrifuge the tube at 230 x g, 20°C for 3 minutes.
6. Remove the supernatant and suspend cells with CellCor™ CD MSC. Count cells and seed to flask at density 4,000-5,000 cells/cm<sup>2</sup> (i.e. AdMSC) with 15 mL of CellCor™ CD MSC.
7. Incubate at 37°C in a humidified atmosphere of 5% CO<sub>2</sub>.

## [Subculture]

1. Remove the cultured medium and wash with 10 mL of PBS.
2. Add 2 mL of Accutase (Life Technologies, Catalog #A1110501) or TrypLE Express (Life Technologies, Catalog #12604013) to the flask and incubate cells for 1-3 minutes at room temperature.  
※ When using detachment reagents such as Trypsin, appropriate inhibitors such as Trypsin inhibitors are required to inactivate enzymes.
3. Observe the cell under the microscope for cell detachment.  
※ Gently tap the flask when cells do not detach.
4. Using CellCor™ CD MSC, collect cells and transfer to the tube.
5. Centrifuge cells at 230 x g, 20°C for 3 minutes.
6. Remove the supernatant and suspend the cell pellets with CellCor™ CD MSC. Count cells and seed to flask at density 4,000-5,000 cells/cm<sup>2</sup> (i.e. AdMSC) with 15 mL of CellCor™ CD MSC.
7. Incubate at 37°C in a humidified atmosphere of 5% CO<sub>2</sub>.
8. Subculture when cells reach up to 75-85% confluence. (Figure 1, usually 3-4 days post seeding)


Figure 1. Cell Confluency (75-85%)

## [Cryopreservation]

1. Remove the supernatant and suspend cell pellets with 1 mL of cryopreservation medium (e.g. Cell-banker2, Masbio, Catalog #11891).
2. Immediately place cryovials in a freezing container and place at -80°C for overnight. After 24 hours, place and store cryovials in a LN<sub>2</sub> tank.

## [Adaptation from medium to CellCor™ CD MSC]

1. After culturing with a prior medium for 24 hours, remove half of the medium and replace the same amount with CellCor™ CD MSC.  
※ i.e. In a T-75 flask, Culture with 15 mL of the prior medium for 24 hours. Remove 7.5 mL of prior medium and add 7.5 mL of CellCor™ CD MSC.
2. When cells reach up to 75-85% confluency, subculture using CellCor™ CD MSC.

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